MI 412 Phosphate Low Range Meter



MARTINI instruments

((

Dear Customer.

Thank you for choosing a Martini product. This manual will provide you with the necessary information for the correct use of the instrument. Please read it carefully before using the meter.

This instrument is in compliance with **C€** directives.

SPECIFICATIONS:	
Range	0.00 to 2.50 mg/L PO ₄
Resolution	0.01 mg/L
Precision	±0.04 mg/L @ 1.00 mg/L
Light Source	Tungsten lamp
Light Detector	Silicon Photocell and 610 nm narrow band interference filter
Method	Adaptation of ascorbic acid method.
Environment	0 to 50°C (32 to 122°F); 100% RH max.
Battery Type	1 x 9 volt
Auto-Shut off	After 10' of non-use
Dimensions	192 x 104 x 52 mm (7.5 x 4.1 x 2")
Weight	380 g

MEASUREMENT PROCEDURE:


1. Turn the meter on by pressing ON/OFF.

10 ml


- 2• When the LCD displays "- -", the meter is ready.
- 3• Fill the cuvet with 10 mL of unreacted sample, up to the mark, and replace the cap.
- 4• Place the cuvet into the holder and ensure that the notch on the cap is positioned securely into the groove.
- 5• Press ZERO and "SIP" will blink on the display.
- 6• After a few seconds the display will show "-0.0-". The meter is now zeroed and ready for measurement.




- 7. Remove the cuvet. Open it.
- 8• Add the content of one packet of MI512-0 reagent.



9• Replace the cap and swirl gently untill the powder is completely dissolved. Reinsert the cuvet into the holder and ensure that the notch on the cap is positioned securely into the groove.





11• The instrument directly displays concentration in mg/L of phosphate.



Interferences

- Iron above 50 mg/L
- Silica above 50 mg/L
- Silicate above 10 mg/L
- Copper above 10 mg/L
- Hydrogen sulfide, arsenate, turbid sample and highly buffered samples also interfere.



GUIDE TO DISPLAY CODES:

This prompt appears for 1 second each time the instrument is turned on.



The dashes "- - -" indicates that the meter is in a ready state and zeroing can be performed.




Sampling In Progress. Flashing "SIP" prompt appears each time the meter is performing a measurement.



"-0.0-", the meter is in a zeroed state and measurement can be performed.



The blinking "BAT" indicates that the battery voltage is getting low and the battery needs to be replaced.



"Conf", the meter has lost its configuration. Contact your dealer or the nearest Martini Customer Service Center.

battery and restart the meter.



ERROR MESSAGES

On zero reading

 Blinking "-0.0-" indicates that the zeroing procedure failed due to a low signalto-noise ratio. In this case press ZERO again.



- "no L", the instrument cannot adjust the light level. Please check that the sample does not contain any debris.
- no L
- "L Lo", there is not enough light to perform a measurement. Please check the preparation of the zero cuvet.



 "L Hi", there is too much light to perform a measurement. Please check the preparation of the zero cuvet.

L Hi

On sample reading

 "-SA-", there is too much light for the sample measurement. Please check if the right sample cuvet is inserted.




• "Inv", the sample and the zero cuvet are inverted.



28-0

 "ZErO", a zero reading was not taken. Follow the instruction in the measurement procedure for zeroing the meter. Under range. A blinking "0.00" indicates that the sample absorbs less light than the zero reference. Check the procedure and make sure you use the same cuvet for reference (zero) and measurement.







ACCESSORIES:

MI0004

MI512-100 Phosphate low range re-

agent (100 tests)

Tissue for wipping cuvets

MI0001 Glass cuvets (2 pcs)
MI0002 Caps for cuvets (2 pcs)
MI0003 Stoppers for cuvets (2 pcs)

(4 pcs)

MI0005 9V battery (1 pc)

BATTERY REPLACEMENT

test.

Battery replacement must only take place in a non-hazardous environment.

Simply rotate the battery cover on the back of the meter.

Detach the battery from the terminals and attach a fresh 9V battery while paying attention to the correct polarity. Insert the battery and replace the cover.



This instrument is warranted against defects in materials and manufacturing for a period of two years from the date of purchase.

If during this period the repair or replacement of parts is required, where the damage is not due to negligence or erroneous operation by the user, please return the meter to either distributor or our office in the original packing and the repair will be free of charge.

Damages due to accidents, misuse, tampering or lack of prescribed maintenance are not covered.

ISTMI412 03/05